


WP3: ctDNA enrichment strategies




This WP is focused on innovative microfluidic-based strategy for processing of patient samples.

Our approach is to combine microfluidic technology with molecular biology to select and enrich ctDNA materials


Task T3.1: Optimization of ctDNA selective capturing on magnetic beads in serum (Curie, M1-M24): (D3.1)

- Streptavidin-biotin system for capture oligo functionalization
- Capture through DNA-DNA interaction
- Improvement of pretreatment throughput and capture efficiency
- Control of hybridization $T^{\circ}\text{C}$



Task T3.2 : Optimization of a ligase-based assay for the selective enrichment of ctDNA targets (Curie, FORTH, M1-M30)

- Oligonucleotides sequence composition and length
- Number of cycles and temperatures
- Use of single-base 3' overhangs
- 5' phosphorylation of oligos
- Introduction of non-complementary tails
- Use of cycling conditions near the oligos T_m .




Task T3.3 : Optimized protocol for ctDNA isolation and enrichment (Curie, FORTH, M1-M36)

- Combination of tasks 3.2 and 3.3 outcomes
- Two strategy direct LCR or DNA release + LCR

Task T3.1: Optimization of ctDNA selective capturing on magnetic beads in serum (Curie, M1-M24): (D3.1)

Lucile Alexandre, Manh Louis Nguyen, Laura Trapiella-Alfonso, Jean-Louis Viovy and Stéphanie Descroix

THE CATCH-U TECHNOLOGY : MICROFLUIDIC FLUIDIZED BED




○ Simplicity of fabrication and operation

Magnetic field gradient

External permanent magnet

MICROFLUIDIC FLUIDIZED BED



- Simplicity of fabrication and operation
- Well controlled particle recirculation and bead density

For Catch-U the beads will be functionalized to extract DNA from patients samples

Task T3.1: ctDNA CAPTURE WITH THE FLUIDIZED BED TECHNOLOGY

Main objectives

1/ Improvement of the fluidized bed technology

✓ Ongoing – 2nd generation

2/ ctDNA capture : two strategies

✓ non specific capture : proof of concept achieved

✓ specific hybridization : ongoing

Fluidized bed 2.0

Main objective: increase the analytical throughput

First generation
H : 50 μm



1 $\mu\text{L/min}$

Generation 2.0
H: 250 μm




1 $\mu\text{L/min}$



5 $\mu\text{L/min}$

Channeling effect

=



Decrease of the ratio of
contact between the liquid
and the solid phase

Necessity to improve the homogeneity
of the bed of beads

- ✓ Bimodal magnetic support
- ✓ Vibration system

Objective 1 : Ongoing work

To improve the fluidized bed throughput:

- Development of the Fluidized bed 2.0
 - ➔ 5 to 20 x Flow rates
 - ➔ Compatible with larger sample volumes
- Homogenization strategies currently investigated : bimodal beads distribution, vibration, chaotic mixing and sequential onf/off state
- Improvement achieved with bimodal sizes of beads and vibration
 - ➔ Capture efficiency - streptavidin/biotin model $> 90\%$


Objective 2 : ctDNA CAPTURE WITH THE FLUIDIZED BED TECHNOLOGY

- ctDNA capture : two strategies
- ✓ non specific capture : proof of concept achieved
- ✓ specific hybridization : ongoing




NON SPECIFIC DNA EXTRACTION


• ChargeSwitch® Technology (Invitrogen)



DNA capture based on electrostatic interactions




NON SPECIFIC EXTRACTION OF FRAGMENTED DNA (100-1000 bp)





All DNA fragments extracted with same efficiency
Potential use for ctDNA

ON-CHIP EXTRACTION EFFICIENCY


200 bp dsDNA



DNA in PBS




DNA in human serum



High capture efficiency of the system despite an important screening effect of serum proteins

DOWNSTREAM ANALYSIS BY DIGITAL PCR



Conclusion and Perspectives

Technology

- Development of fluidized 2.0 able to accommodate large volume
- Development of different strategies to improve the bed homogeneity

DNA capture

- Demonstration of non specific DNA capture on chip with artificial samples and patients samples
- Preliminary results on specific DNA capture

DNA capture combined with LCR

- Preliminary results on DNA release
- Strategy for LCR on chip – Exchange of PhD in July with FORTH